



Cannabis Manufacturing & Testing



Supplies • Chemicals • Safety • Equipment

...The Gem of Lab Distributors

734-405-6620 • Customerservice@jadesci.com

Partnering with Precision: Jade Scientific, Inc.

As the cannabis industry evolves from emerging to established, the demand for precision, compliance, and reproducibility in testing and production has never been higher. At the intersection of science and innovation stands Jade Scientific, Inc.—a trusted partner supplying the tools cannabis manufacturers and laboratories need to operate efficiently, accurately, and within regulatory frameworks.

Your Trusted Laboratory Solutions Provider:

Founded on a passion for science and service, Jade Scientific, Inc. has become a go-to source for laboratories across the life sciences, diagnostics, and cannabis sectors. Whether you're outfitting a new lab or upgrading your current capabilities, Jade offers a broad portfolio of essential lab equipment, consumables, and instrumentation tailored for cannabis testing workflows.

Essential Solutions for Cannabis Testing and Manufacturing:

Analytical Instruments:

- HPLC and GC systems for cannabinoid and terpene quantification
- LC-MS and GC-MS platforms for pesticide, solvent, and heavy metal analysis

Sample Preparation:

- QuEChERS kits, SPE cartridges, homogenizers, and centrifuges for efficient sample prep
- Precision balances, filtration systems, and mixers to ensure batch consistency

Processing Equipment:

- Washdown scales, ultra-low freezers, freeze dryers, drying ovens, winterization kits & components, distillation systems & glassware
- Mobile shelving & drying racks, moisture analyzers, water activity meters, bulk food grade storage solutions

Certified Reference Standards:

- CRMs for cannabinoids, terpenes, residual solvents, and pesticides
- Designed for ISO/IEC 17025 and GMP-compliant operations

Environmental and Safety Products:

- Fume hoods, PPE, cleanroom supplies, and surface disinfectants
- Maintain safety and reduce contamination risks in lab and production environments



Jade Scientific is proud to operate under an ISO 9001:2015 compliant quality management system, reflecting our commitment to consistent quality, customer satisfaction, and continuous improvement. This internationally recognized standard ensures that every aspect of our operations—from sourcing and inventory to order fulfillment and customer support—meets rigorous quality benchmarks. For our customers in the cannabis industry, it means reliable service, dependable products, and a partner you can trust.

Supplying the Industry's Most Trusted Brands

At Jade Scientific, we pride ourselves on offering products from the most trusted and respected brands in the laboratory industry. We carefully curate our catalog to include manufacturers known for quality, reliability, and regulatory compliance—ensuring our customers have access to equipment and consumables they can depend on. From industry leading chromatography columns and certified reference standards to precision balances, glassware, and safety supplies, our partnerships with top-tier brands allow us to support cannabis testing and manufacturing operations with confidence and consistency. Whether you're outfitting a new lab or upgrading your current workflows, you can trust Jade Scientific to deliver the best names in the business.

We offer products from the following manufacturers and many others:



Safety is a top priority in every stage of cannabis manufacturing and testing. Jade Scientific, Inc. provides essential safety products including PPE, spill control materials, and chemical storage solutions to help protect personnel and maintain a compliant lab environment. Our safety offerings support both routine operations and emergency preparedness. With trusted brands and reliable inventory, we help ensure your facility meets industry standards and operates with confidence.

Epic® Tian's International, Inc. Protection Coverall

Latex free / low lint coverall for use in some critical / clean environments. Nonwoven fabric is cloth-like and cool to wear. Features include an attached hood, attached boots, zipper front entry and elastic wrists.



Part No.	Description	Unit
TIA-200881-L	Coveralls, White Polypro, H&B, EW, Large	25
TIA-200881-2XL	Coveralls, White Polypro, H&B, EW, 2X Large	25
TIA-200881-3XL	Coveralls, White Polypro, H&B, EW, 3X Large	25
TIA-200881-4XL	Coveralls, White Polypro, H&B, EW, 4X Large	25
TIA-200881-5XL	Coveralls, White Polypro, H&B, EW, 5X Large	25

Epic® Tian's International SafeTrack Shoe Covers

Latex free / lint free cleanroom boot cover with seamless bottom. Unique laminated fabric offers unbelievable traction and skid resistance. This durable fabric is also impervious to liquids and suitable for use in clean environments. Elasticized top opening.



Part No.	Description	Unit
TIA-537782-L	Shoe Cover, SAFETRACK High Traction, Blue, L	200
TIA-200881-2XL	Shoe Cover, SAFETRACK High Traction, Blue, 2XL	200

FyterTech Nonwovens Poly-Drum Hazmat Spill Kits

- 20-, 30-, 55- and 95-gallon kits
- Absorb up to 16-, 24-, 40-, or 82-gallons
- Aggressive and non-aggressive spills



Part No.	Description	Unit
NPS-250020	20 Gallon Hazmat Spill Kit (Absorbs Up To 16 Gal)	1
NPS-250030	30 Gallon Hazmat Spill Kit (Absorbs Up To 24 Gal)	1
NPS-250055	55 Gallon Hazmat Spill Kit (Absorbs Up To 40 Gal)	1
NPS-250095	95 Gallon Hazmat Spill Kit (Absorbs Up To 82 Gal)	1

Hourglass International Nitrile Exam Gloves

HandPRO® FreeStyle1100™ Nitrile Exam Gloves with Low Dermatitis Potential are highly flexible and conform to the hand for exceptional freedom of movement. Made from the patented LOW DERMATM accelerator-free nitrile, this powder-free glove has fewer chemicals and virtually eliminates the risk of Type I latex and Type IV chemical allergic reactions and contact dermatitis. The unique nitrile compound sheds water better for a more secure grip – even in wet environments.



Part No.	Description	Unit
HGI-1100	FreeStyle1100™ Nitrile Exam Gloves XS	1000
HGI-1101	FreeStyle1100™ Nitrile Exam Gloves S	1000
HGI-1102	FreeStyle1100™ Nitrile Exam Gloves M	1000
HGI-1103	FreeStyle1100™ Nitrile Exam Gloves L	1000
HGI-1104	FreeStyle1100™ Nitrile Exam Gloves XL	1000

Hourglass International Nitrile Exam Gloves Extended Cuff

HandPRO® Fortis500™ Extended Cuff Powder-Free Nitrile Exam Gloves are created from the latest patented LOW DERMA™ technology* accelerator-free nitrile. The glove is thicker than typical exam gloves for increased puncture and abrasion resistance, and features an extended cuff for extra wrist and forearm protection in demanding, high-risk environments.



Part No.	Description	Unit
HGI-500	Fortis500™ Extended Cuff Nitrile Exam Gloves S	1000
HGI-501	Fortis500™ Extended Cuff Nitrile Exam Gloves M	1000
HGI-502	Fortis500™ Extended Cuff Nitrile Exam Gloves L	1000
HGI-503	Fortis500™ Extended Cuff Nitrile Exam Gloves XL	1000
HGI-504	Fortis500™ Extended Cuff Nitrile Exam Gloves XXL	1000

Find a larger variety of gloves and brands online at jadesci.com

High-Tech Conversions, Inc Cleanroom Entrance Mats

Tacky Traxx™ Mats are made from low density polyethylene sheets with an adhesive for trapping dirt and particles from foot or wheel traffic. These mats will ensure cleanliness upon entry into controlled environment.

- Effectively captures dirt and dust from both foot traffic and equipment wheels prior to entering cleanroom.
- Removal of sheet peels off smoothly to reveal a new clean surface.
- Each sheet is numbered allowing user to track remaining sheets with ease.
- Smooth and effective sheets is not damaged by shoes, booties or cart wheels.



Part No.	Description	Unit
HTC-TT-1836-34W	White, 18" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-1836-34B	Blue, 18" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-1836-34G	Grey, 18" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-1845-34W	White, 18" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-1845-34B	Blue, 18" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-1845-34G	Grey, 18" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2436-34W	White, 24" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-2436-34B	Blue, 24" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-2436-34G	Grey, 24" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-2445-34W	White, 24" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2445-34B	Blue, 24" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2445-34G	Grey, 24" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2645-34W	White, 26" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2645-34B	Blue, 26" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-2645-34G	Grey, 26" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-3636-34W	White, 36" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-3636-34B	Blue, 36" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-3636-34G	Grey, 36" X 36" Sticky Mat	4 X 30 Sheets
HTC-TT-3645-34W	White, 36" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-3636-34B	Blue, 36" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-3636-34G	Grey, 36" X 45" Sticky Mat	4 X 30 Sheets
HTC-TT-3660-34W	White, 36" X 60" Sticky Mat	4 X 30 Sheets
HTC-TT-3660-34B	Blue, 36" X 60" Sticky Mat	4 X 30 Sheets
HTC-TT-3660-34G	Grey, 36" X 60" Sticky Mat	4 X 30 Sheets

Epic® Tian's International, Inc. 3 Pocket Lab Coats

Latex free / low lint lab coat for use in controlled environments and many general application areas. Breathable cloth-like fabric is comfortable and cool. Features include plastic snaps, folded collar, knit wrists and 3 pockets.



Part No.	Description	Unit
TIA-845884-M	Basic Protection CleanRoom Lab Coat, M	1
TIA-845884-L	Basic Protection CleanRoom Lab Coat, L	1
TIA-845884-XL	Basic Protection CleanRoom Lab Coat, XL	1
TIA-845884-2XL	Basic Protection CleanRoom Lab Coat, 2XL	1

Heathrow Scientific CRYO-PROTECTION® Face Shield

- The injection molded polycarbonate clear shield is shaped to hug your chest and side of your face when in use, and pivots up and out of the way when your task is done.
- The padded head gear has a ratchet band in the back and a pin-and-lock system over the top for a 2-way adjustable fit.
- A lightweight design meets ANSI Z87.1 requirements, is CE Cat. II, and is impact resistant.
- Rigorous testing has confirmed that the shield retains its optical clarity and impact resistance after being splashed with liquid cryogens.



Part No.	Description	Unit
HRS-FSH1001HG	CRYO-PROTECTION® Face Shield Headgear	1
HRS-FSH1001WIN	CRYO-PROTECTION® Face Shield Window	1
HRS-FSH1001	CRYO-PROTECTION® Face Shield Complete	1

Epic® Tian's International, Inc. Polypropylene Beard Covers

Latex free, low lint, 100% machine made cover expands from .75 x 5" packed size and opens up to a full sized beard cover. Machine welded seam eliminates all needle holes, thread and stitching. Double elastic faceband for a secure and comfortable fit. Lightweight and cool/breathable non-woven fabric.



Part No.	Description	Unit
TIA-91588	Beard Cover	1000

Epic® Tian's International, Inc. Cleanroom Pleated Bouffant Cap

Latex free, low lint full coverage cleanroom bouffant style cap. 100% machine made cap expands from 1"x7" packed size and opens up to a full sized cap. Machine welded seam eliminates all needle holes, thread and stitching. Double elastic headband for a secure and comfortable fit. Lightweight and cool/breathable non-woven fabric



Part No.	Description	Unit
TIA-31588-21	CRYO-PROTECTION® Face Shield Headgear	1000
TIA-31588-24	CRYO-PROTECTION® Face Shield Window	1000

Dynalab Corp Glass Disposal Box

- Floor model glass disposal box with handles is safer in use and disposal
- Glass Disposal Boxes feature colorful graphics for quick identification
- Assembly instructions on each box
- Tough 2 mil polyethylene liner
- Folding design for shipping and storage
- Closure assures a safe, hygienic disposal
- Glass Disposal Boxes come in a range of convenient sizes
- Replacement PE bags available for Floor Model Glass Disposal Boxes



Part No.	Description	Unit
DYN-797015	Glass Disposal Bench Box	6
DYN-797035	Glass Disposal Box Intermediate	6

Eisco Scientific Safety Goggles

Designed to protect your eyes from chemical splashes, high speed particles, dust and UV wavelengths (360-380). Lightweight, durable polycarbonate lens. The soft, flexible vinyl frame and pliable nose bridge provides a comfortable fit. Designed to be able to fit over prescription eyeglasses. Goggles feature a freely adjustable, soft elastic head strap.



Part No.	Description	Unit
EIS-G103	Basic Green Goggle	1
EIS-G103CLR	Basic Clear Goggle	1
EIS-G104	Vented Green Goggle	1

Eisco Scientific Safety Glasses

Durable and lightweight, impact resistant polycarbonate lens with UV protection. Conforms to ANSI Z87.1 2010 sideshield protection. Side vents allow for airflow and prevent lenses from fogging. Designed to be able to fit over prescription eyeglasses. Size fits most faces comfortably and securely.



Part No.	Description	Unit
EIS-CH0619A	Basic Green Goggle	1

Bel-Art Emergency Eye Wash Stations

- Three stations are available to suit all workstations: a single or a double station with 1000ml (32 oz) eyewash bottle(s) or a single station with a 500ml (16 oz) bottle
- Shipped with empty bottles; can be filled with your favorite Sterile Eye Wash Solution Refills
- Quick instruction diagrams on the bottle and station confirm how to use the eyewash bottle
- Includes a mirror and multi-language directions on the station (English, French, German and Spanish)
- Single station dimensions: 44.5W x 38cmH (17½ x 15"); double station dimensions: 45.7W x 34.9cmH (18 x 13¾")
- Hardware not included



Part No.	Description	Unit
BEL-24865-0000	Eyewash station with 16 oz. bottle	1
BEL-24868-0000	Double Eyewash station with two 16 oz. bottles	1
BEL-24845-0000	16 oz. LDPE Eye Wash Bottle	1

3M Half Facepiece Respirator

- Economical and low maintenance
- Lightweight elastomeric material for comfort
- Low-profile design minimizes interference with field of view



Part No.	Description	Unit
3MP-6200	3M Half Facepiece Reusable Respirator, Medium	1

Thermo Scientific Barnstead™ Pacific™ RO Water Purification System

Improve functional and economical production of high-purity water with this modularly designed laboratory reverse osmosis system.

Part No.	Description	Unit
TFB-50132385	Pacific RO 3	1
TFB-50132386	Pacific RO 7	1
TFB-50132387	Pacific RO 12	1
TFB-50132388	Pacific RO 20	1
TFB-50132389	Pacific RO 40	1



ELGA Lab Water PURELAB® Flex

Bringing you more of the flex you demand with a significantly reduced environmental impact. Providing perfect water purity for applications requiring ultrapure type I (18.2 MΩ.cm) water, the PURELAB flex allows you to focus on what you do best.



Part No.	Description	Unit
ELG-SK-PF2XXXXM1-Kit	PURELAB Flex 2, Water Purification System, with starter kit.	1

Millipore Sigma Milli-Q® IQ 7003/05/10/15 Ultrapure & Pure Lab Water Purification System

This fully-integrated water purification system is your ultimate lab water solution providing:

- Superior and consistent quality Type 1 and Type 2 water directly from tap water
- No maintenance or hazardous chemicals with Elix® EDI technology
- Ergonomic Q-POD® and E-POD® dispensers
- Greater confidence in lab productivity with MyMilli-Q™ digital services
- A reduced environmental footprint with the mercury-free UV lamps and less electricity, materials and packaging



Part No.	Description	Unit
EMD-IPAKKITA1	Milli-Q® IQ/EQ 7000 Purification kit	1
EMD-ZIQ7000T0C	Milli-Q® IQ/EQ 7000 Purification System	1

Ohaus DEFENDER™ Series Bench Scale

The Defender Series features a simple, yet rugged tubular frame base design and indicators with tactile keys, multi-color backlit LCD display, AC or battery operation, and flexible mounting capabilities.



Part No.	Description	Unit
OHA-30461648	DEFENDER® 5000 - D52 Multifunctional Bench Scale	1
OHA-30419640	DEFENDER® 5000 Base	1

Benchmark Scientific BeadBlaster™ 96 Ball Mill Homogenizer

The BeadBlaster™ 96 is an extremely versatile bead mill homogenizer that has applications in a variety of areas, including biological research, environmental testing, and industrial settings. Adapters are available for microplates, microtubes, and 50ml tubes.



Part No.	Description	Unit
BEN-IPD9600	BeadBlaster 96 Ball Mill Homogenizer, 120V US Plug	1
BEN-IPD9600-E	BeadBlaster 96 Ball Mill Homogenizer, 230V EU Plug	1
BEN-IPD9600-96	Adapter for 96 well plate	2
BEN-IPD9600-96A2	Tube Rack, 20 x 2ml, Pair, for BeadBlaster 96. Requires IPD9600-96	2

Benchmark Scientific Overhead Stirrers

The BeadBlaster™ 96 is an extremely versatile bead mill homogenizer that has applications in a variety of areas, including biological research, environmental testing, and industrial settings. Adapters are available for microplates, microtubes, and 50ml tubes. Grinding jars and balls in a variety of materials are available.



Part No.	Description	Unit
BEN-IPS2050-20	OS20L Overhead Stirrer, 20L includes stand, and 4 arm propeller (IPS2050-P-S4), 100 to 240V with USA Plug	1
BEN-IPS2050-20-E	OS20L Overhead Stirrer, 20L includes stand, and 4 arm propeller (IPS2050-P-S4), 100 to 240V with EU Plug	1

Benchmark Scientific V2™ Vortex Mixer

- Reversible top platform, cup or flat-dimpled surface
- Powerful motor for INSTANT vortexing
- Continuous or “touch” operation
- Q-Drive™ dynamic balancing system
- Variable speed from 200 to 3200 rpm
- Wide variety of accessories, for tubes up to 50ml.



Part No.	Description	Unit
BEN-BV1003	BenchMixer V2™ Vortex Mixer with flip top cup head and new counter balance, 115V	1
BEN-BV1003-E	BenchMixer V2™ Vortex Mixer with flip top cup head and new counter balance, 230V	1

Heidolph Hei-VAP Core Rotary Evaporator

The BeadBlaster™ 96 is an extremely versatile bead mill homogenizer that has applications in a variety of areas, including biological research, environmental testing, and industrial settings. Adapters are available for microplates, microtubes, and 50ml tubes.



Part No.	Description	Unit
HEI-036040065	Hei-VAP Core ML/ G1B Coated Rotary Evaporator	1
HEI-036040066	Hei-VAP Core ML/ G3B Coated Rotary Evaporator	1
HEI-036040068	Hei-VAP Core ML/ G5B Coated Rotary Evaporator	1
HEI-036040069	Hei-VAP Core ML/ G6B Coated Rotary Evaporator	1

Labconco Scroll Vacuum Pump

Scroll Vacuum Pumps are low maintenance, corrosion resistant pumps that require no oil. It is recommended that the scrolls be changed every 40,000 hours of use. Eliminating the oil eliminates hydrocarbon contamination. Scroll pumps handle water vapor better than most types of pumps and make significantly less noise during operation.



Part No.	Description	Unit
LCC-7587000	Scroll Vacuum Pump, 190 liters/min, 115V, 50/60 Hz	1
LCC-7587100	Scroll Vacuum Pump, 103 liters/min, 115/230V, 50/60Hz	1
LCC-7587300	Silencer for Scroll Vacuum Pump	1

Branson Brasonic® Ultrasonic Baths

- Easily accessible panel located above and behind the tank for safety and reliability, to prevent spills on the control panel, back lit on the digital models
- Built-in drains on the 1.5 gallon units and larger, .5 and .75 gallon have a convenient built-in pour drain
- Sweep frequency
- 40KHz industrial transducers



Part No.	Description	Unit
BRA-CPX-952-119R	Branson Digital ultrasonic, 1.9 L (0.5 gal.), 120V	1
BRA-CPX-952-219R	Branson Digital ultrasonic, 2.8 L (0.75 gal.), 120V	1
BRA-CPX-952-319R	Branson Ultrasonic Bath, 5.7L (1.5 gallons), 120V	1
BRA-CPX2800H	Branson Ultrasonic Cleaning Bath, 0.75 gallon	1
BRA-CPX-952-216R	Ultrasonic Cleaner with Mechanical Timer and Tank Capacity of 3/4 Gallon (2.8L) 120 V, 50/60 Hz	1

PolyScience Digital General Purpose Water Baths

PolyScience water baths are the first baths on the market to utilize a full color TFT display. Our clear hinged gable cover continues to be a standard feature on all of our water baths and the intuitive plain-language user interface offers five selectable languages.



Part No.	Description	Unit
PLY-WBE02A11B	2L Digital Water Bath	1
PLY-WBE05A11B	5L Digital Water Bath	1
PLY-WBE10A11B	10L Digital Water Bath	1
PLY-WBE20A11B	20L Digital Water Bath	1

Julabo CORIO laboratory circulators

JULABO circulators are used successfully around the world. Whether in research, material testing or in production – the well proven and reliable technology is valued by users in all industries. By focusing on user requirements JULABO has set the benchmark for temperature control technology for over 4 decades.



Part No.	Description	Unit
JUL-9012413.02	CORIO CD-B13 115V/60Hz	1

Larger ovens and options available at Jadesci.com

Ohaus

Analytical Balance PR224/E

Offering accuracy and repeatability in essential weighing applications in laboratory, industrial and education settings, PR Analytical and Precision Balances deliver competitive performance at an economical price. Featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation.



Part No.	Description	Unit
OHA-30430062	Analytical Balance PR224/E	1

Ohaus

Guardian Hotplates Stirrers

- Temperature range up to +80°C
- 8 model sizes from 50 up to 750 liters
- 2 variants: SingleDisplay and TwinDisplay
- Temperature ramps programmable via AtmoCONTROL software (models with TwinDISPLAY only)



Part No.	Description	Unit
OHA-30500530	Guardian 5000 HP 7x7 Ceramic	1
OHA-30500550	Guardian 5000 HPS 10x10 Ceramic	1
OHA-30500560	Guardian 5000 HPS Round Aluminum	1
OHA-30500540	Guardian 5000 ST 7x7 Ceramic	1
OHA-30500610	Guardian 7000 HPS 10x10 Ceramic	1
OHA-30500600	Guardian 7000 HPS 7x7 Ceramic	1
OHA-30500620	Guardian 7000 HPS Round Aluminum	1
OHA-30500520	Guardian 5000, 7x7 Ceramic, (120V)	1

Benchmark Scientific

Benchtop Multitube Vortex

- Up to 2500rpm, a true vortex in every tube
- Assortment of accessories for a variety of tubes, flasks & plates
- Pulsing program for interval mixing
- Fits almost any tube rack



Part No.	Description	Unit
BEN-BV1010	BenchMixer™ XL Multi-Tube Vortexer, includes 50x12mm tube rack, 115V	1
BEN-H1010-E	BenchMixer™ XL Multi-Tube Vortexer, includes 50x12mm tube rack, 230V	1

Benchmark Scientific

MyTemp H Mini Digital Incubators

- Digital temperature control
- Personal sized, 13 x 14.5" footprint (20 L)
- Accepts bottles and flasks up to 2 L
- Two models: Heat Only or Heat/Cool



Part No.	Description	Unit
BEN-H2200-H	MyTemp™ Mini Digital incubator	1
BEN-H2200-HC	MyTemp™ Mini Digital Incubator, Heating/Cooling 115V	1
BEN-H2200-HC-E	MyTemp™ Mini Digital Incubator, Heating/Cooling 230V	1
BEN-H2200-H-E	MyTemp™ Mini Digital Incubator, Heating only 230V	1

Ohaus

Moisture Analyzers

Features a redesigned housing with stronger material, innovative carbon fiber heating technology and an upgraded load cell for enhanced durability and performance. The easy-to-use MB32 is the perfect choice for routine moisture testing.



Part No.	Description	Unit
OHA-83999801	Moisture analyzer, Model MB-23	1
OHA-30251683	Moisture Analyzer, MB27 AM	1
OHA-30241165	Moisture Analyzer, MB120 AM	1

Benchmark Scientific

Hermle High Capacity Centrifuges

- High capacity swing out rotor, 68 x 15ml or 28 x 50ml
- Ideal for tissue culture applications
- Advanced programming with automatic rotor recognition
- Rotor exchange is simple & requires less than 1 minute
- Available with refrigeration



Part No.	Description	Unit
BEN-Z446	Z446 Centrifuge, 120V	1
BEN-Z446-E	Z446 Centrifuge, 230V	1
BEN-Z446-K-E	Z446-K Refrigerated Centrifuge, 230V	1

Benchmark

Hermle Ultra High Capacity Centrifuges

- Ultra High Capacity (4x1000ml, 98x15ml or 40x50ml)
- Accepts 500ml or 750ml blood bags
- Two models:
 - Floor-standing, with storage drawer
 - Underbench (27.5" height)



Part No.	Description	Unit
BEN-Z496-K	Z496-K High Volume Floor Standing centrifuge	1
BEN-Z496-K-UC	Z496-K High Volume Underbench centrifuge	1

Benchmark

Digital Hotplate Magnetic Stirrer

- Large, ceramic work surface, 7 x 7 in
- Safety LED indicates hot surface
- 5L stir capacity
- Control actual sample temperature (with optional probe)



Part No.	Description	Unit
BEN-H3770-H	Digital Hotplate, 115V	1
BEN-H3770-H-E	Digital Hotplate, 230V	1
BEN-H3770-HS	Digital Hotplate Stirrer, 115 V	1
BEN-H3770-HS-E	Digital Hotplate Magnetic Stirrer, 230V	1
BEN-H3770-S	Magnetic Stirrer, 115V	1
BEN-H3770-SE	Magnetic Stirrer, 230V	1

Quincy Labs

Gravity Convection Ovens

- Offered in four different sizes to fit your needs.
- Gravity convection (GC) or Air Forced (AF).
- Hard scratch-resistant hammer finish exterior.
- Corrosion-resistant aluminized steel Interior.
- One inch of high-density mineral wool insulation.
- Temperature is controlled by a hydraulic thermostat.
- UL / CUL listed.

Part No.	Description	Unit
QLB-20GC	Lab Oven, 1.27 cubic feet, Gravity Convection	1
QLB-30GC	Gravity Convection Oven 2cu ft,max 232C/450F	1
QLB-30GCE	Gravity Convection Oven 2cu ft,max 232C/450F	1



Labconco

Protector XL Floor-Mounted Fume Hoods

Protector XL Floor-Mounted (Walk-In) Laboratory Hoods have panel liners that may be constructed of a variety of corrosion-resistant materials and supplied in three depths to meet individual installation and application requirements. Suitable for general chemistry applications and accommodate oversized apparatus. Widths of up to 16 feet are offered.

- Dry powder, glacier white epoxy-coated steel exterior
- Chemical-resistant molded composite panel liner and adjustable upper and lower baffles with flame spread less than 25 per ASTM E-84
- Removable front and side panels and exterior and interior service access panels for access to plumbing and electrical wiring
- Pre-wired instant-start T8 fluorescent lighting, light and blower switches for 115 volt, 60 Hz operation
- Tempered safety glass vertical-rising sashes framed with PVC (4', 5' and 6' models) or epoxy-coated aluminum (8' models)
- Tempered safety glass horizontal-sliding sashes with epoxy-coated aluminum frames and sash tracks (8', 10', 12', 16' models)



Part No.	Description	Unit
LCC-112836000	8' Protector XL Floor-Mounted Hood with horizontal-sliding sashes, 55.7" exterior depth, 100-115V, 50/60 Hz	1
LCC-112836002	8' Protector XL Floor-Mounted Hood with horizontal-sliding sashes, 55.7" exterior depth, with 2 service fixtures and 2 GFCI duplexes, 100-115V, 50/60 Hz	1
LCC-112836020	8' Protector XL Floor-Mounted Hood with horizontal-sliding sashes, 55.7" exterior depth, 208-230V, 50/60 Hz	1
LCC-112836021	8' Protector XL Floor-Mounted Hood with horizontal-sliding sashes, 55.7" exterior depth, with 2 service fixtures, 208-230V, 50/60 Hz	1
LCC-112837000	8' Protector XL Floor-Mounted Hood with vertical-rising sashes, 55.7" exterior depth, 100-115V, 50/60 Hz	1
LCC-112837002	8' Protector XL Floor-Mounted Hood with vertical-rising sashes, 55.7" exterior depth, with 2 service fixtures and 2 GFCI duplexes, 100-115V, 50/60 Hz	1
LCC-112837020	8' Protector XL Floor-Mounted Hood with vertical-rising sashes, 55.7" exterior depth, 208-230V, 50/60 Hz	1
LCC-112837021	8' Protector XL Floor-Mounted Hood with vertical-rising sashes, 55.7" exterior depth, with 2 service fixtures, 208-230V, 50/60 Hz	1

Benchmark

MyFuge™ Mini Centrifuge

- COMBI-Rotor for tubes and strips
- Twice the capacity of most mini centrifuges
- Near silent operation
- Starts & stops with closing/opening of the lid



Part No.	Description	Unit
BEN-C1008	MyFuge™ Mini Centrifuge, 2 rotors, 100-240V	1
BEN-C1012	MyFuge™12 Mini Centrifuge 100-240V	1

Thermo Fisher Scientific AccuTemp Vacuum Drying Ovens

Dry up to six times faster than when using conventional units. These ovens feature excellent heat transfer and fast heating up time, minimizing operation time.



Part No.	Description	Unit
TFB-51023262	Vacutherm VT6130M Vacuum Oven, jacket heated, 128L Volume, 208V, 60 Hz	1
TFB-51023263	Vacutherm VT6060P Vacuum Oven, shelf heated, 53L volume, 120V, 60 Hz	1
TFB-51023264	Vacutherm VT6130P Vacuum Oven, shelf heated, 128L Volume, 208V, 60 Hz	1

United Scientific Supplies Heating Mantles, Combo

Units feature space-saving built-in electronic power regulator knobs as well as polyurethane-coated aluminum exterior and high-quality knitted mantles with glass fiber underlay for uniform distribution of heat. Maximum operating temperature is 400°C; mantles can be used in a wide range of ambient temperatures, from about 45°F to 100°F, and relative humidity up to 80%. Mantles are designed to securely hold flasks so that the top half of the flask is easily visible during operation.



Part No.	Description	Unit
USS-HM3P2000-100	Heating Mantles, Combo, 100MI	1
USS-HM3P2000-250	Heating Mantles, Combo, 250MI	1
USS-HM3P2000-500	Heating Mantles, Combo, 500MI	1
USS-HM6P2000-250	Heating Mantles, Combo, 250MI	1

Benchmark Scientific D1000 Homogenizer

- Homogenize most samples in under 30 seconds
- For volumes from 100µl to 250ml
- Supplied with both 5mm & 7mm generator probes



Part No.	Description	Unit
BEN-D1000	D1000 Homogenizer, includes 5mm and 7mm generators, 115V	1
BEN-D1000-E	D1000 Homogenizer, includes 5mm and 7mm generators, 230V	1

BrandTech® Transferpette Electronic Single Channel Pipettes with Charger

The BRAND Transferpette® electronic is an ergonomic, intuitive electronic pipette for use in most laboratory applications. It is available in five volume ranges for reliable dispensing of volumes from 0.5µl to 5ml. The soft-touch pipetting key, adjustable finger-rest, light weight and innovative tip cone minimize static strain and operational forces, reducing the risk of repetitive stress injuries (RSI).



Part No.	Description	Unit
BRD-705319	Transferpette electronic single channel 0.5-10uL	1
BRD-705320	Transferpette electronic single channel 2-20uL	1
BRD-705323	Transferpette electronic single channel, 20-200uL	1
BRD-705326	Transferpette electronic single channel 100-1000uL	1
BRD-705327	Transferpette electronic single channel 500-5000uL	1

Haier -86°C Ultra-Low Temperature Freezers

- Microprocessor-controlled system designed for controllable range of -40°C to -86°C for cabinet space with 1°C increment
- Large LED display for cabinet temperature, set temperature, ambient temperature, and input voltage
- Settable high temperature and low temperature alarms
- Automatic clean-filter alarm and sensor error alert
- Adjustable storage shelf height
- Can be connected with network, and remote alarm systems



Part No.	Description	Unit
HAI-DW-86L100J	Ultralow Freezer; -86°C, Ultra-Low, 3.5 cu. ft., Undercounter, Haier Biomedical, 120V	1
HAI-DW-86L728BPST	Ultra-Low Temperature -86C Freezer, Energy Star 728L (26cf), 115V/60Hz	1
HAI-DW-86L728J	Ultra-Low Temperature -86C Freezer, Energy Star 728L (26cf), 115V/60Hz	1
HAI-DW-86L728ST	Ultralow Dual Autocascade -86C Freezer, Energy Star 728L (26cf), 115V/60Hz	1
HAI-DW-86L729BPT	Ultralow Freq. Conversion (variable speed) -86C Freezer, Energy Star 729L (26cf) 110-230V/60Hz	1
HAI-DW-86L829BPT	Ultralow Freq. Conversion (variable speed) -86C Freezer, Energy Star 829L (29cf) 110-230V/60Hz	1
HAI-DW-86L959BPT	Ultralow -86C Freq. Conversion (variable speed) freezer, Energy Star 959L (34f), 208-230V/60Hz	1

Heidolph Tuttnauer Vertical Autoclaves

- No utility connections needed (models 2840 ELV-WR and 3870ELV-WR) just plug in and add water
- Two integrated temp. probes to ensure the liquid is completely sterilized
- Four validated cycles (glass, plastic, liquid A and liquid B-waste) and two liquid cooling cycles included with 5075ELVC
- Stores up to 200 past cycles you can open on your PC
- Create up to 16 custom cycles
- Two baskets included
- Optional biohazard filter (BH) prevents contaminated air from exhausting into the lab



Part No.	Description	Unit
HEI-023210990	Heidolph Tuttnauer Lab Line Vertical Autoclave 5075ELVC	1
HEI-023210744	Heidolph Tuttnauer Lab Line Vertical Autoclave 2840ELV	1
HEI-023210745	Heidolph Tuttnauer Lab Line Vertical Autoclave 3870ELV	1

Labconco Fume Hood

All fume hoods include a 3/16" (0.5 cm) thick tempered safety glass counterbalanced sash and a 10" (25.4 cm) diameter duct collar. Hoods are constructed of glacier white epoxy-coated 16 gauge steel and include a removable front panel and a 4" (10.2 cm) L x 4 1/2" (11.4 cm) W utility pass-through port with molded neoprene cover.



Part No.	Description	Unit
LCC-2247300	Basic 47 Hood with blower, 115V, 60 Hz	1
LCC-2247301	Basic 47 Hood with blower, 220V, 50 Hz	1
LCC-2247400	Basic 47 Hood without blower, 115V, 60 Hz	1
LCC-2247401	Basic 47 Hood without blower, 220V, 50 Hz	1
LCC-2247500	Basic 47 Hood with EP blower	1
LCC-2247501	Basic 47 Hood with EP blower, 220V, 50 Hz	1
LCC-2246300	Basic 70 Hood with blower, 115V, 60 Hz	1
LCC-2246301	Basic 70 Hood with blower, 220V, 50 Hz	1
LCC-2246400	Basic 70 Hood without blower, 115V, 60 Hz	1
LCC-2246401	Basic 70 Hood without blower, 220V, 50 Hz	1
LCC-2246500	Basic 70 Hood with EP blower, 115V, 60 Hz	1
LCC-2246501	Basic 70 Hood with EP blower, 220V, 50 Hz	1

Ohaus Rocking Shaker

OHAUS Rocking and Waving Shakers are innovatively designed for use in a wide range of lab applications. Rocking Shakers provide controlled seesaw-like motion, while Waving Shakers offer a smooth, low-foaming, 3D wave motion. Available in microprocessor-controlled digital models and economical non-digital models. All models include a non-slip mat, and most models are available with a stacking tray option to easily double the capacity.



Part No.	Description	Unit
OHA-30391973	Shaker, Waving, SHWV02AL, US	1
OHA-30391994	Shaker, Rocking, SHRK04DG, US	1

JustRight Acid Storage Cabinets

- Reinforced double wall construction with 1-1/2 in (3.8 cm) air space on sides, top, bottom, back and doors
- Constructed of 18-gauge steel
- Lead-free, high-gloss, baked-on powder finish.
- Continuous piano hinges
- Ground connection
- Self-latching doors



Part No.	Description	Unit
EAG-CRA32	Corrosive Storage Cabinet, 30 gallon, Manual Door, Blue	1

More sizes and options available at Jadesci.com

Labconco Fiberglass Corrosive-Resistant Blower

Fiberglass Blowers are available that deliver as little as 300 CFM at .12" static pressure or up to 2050 CFM at 4.0" static pressure. Fiberglass Blowers are not recommended for use with Perchloric Acid.



Part No.	Description	Unit
OHA-30391973	Shaker, Waving, SHWV02AL, US	1
OHA-30391994	Shaker, Rocking, SHRK04DG, US	1

Eisco

Distillation Kit Sets

- Sets are made of high quality Borosilicate glass.



Part No.	Description	Unit
EIS-CH1026	Distillation Set, Simple, 24/29 Joints, 4pc	1
EIS-CH0879	Distillation Kit, Micro Glass, 25ml, 9 Piece Set	1

More available at Jadesci.com

PolyScience

Refrigerated Circulator 7 Liter

- Intuitive 3.75" (9.5 cm) display with touch-pad control
- 4 languages: French, German, Spanish, English
- Swivel 180™ Rotating Controller
- On-board connectivity: USB-A & B, Ethernet, RS232/addressable RS485, and external temperature probe
- On-screen prompts
- Automatic performance optimization and specific heat tuning



Part No.	Description	Unit
PLY-AD07R-20-A11B	7L Space-Saving Refrigerated Circulator, Advanced Digital (-20° to 200°C), 120V, 60Hz	1

Eisco Scientific

Goggle Sanitizer Cabinet

Place up to 36 pairs of goggles into this goggle cleaning cabinet. (Goggles are NOT included.) Two safety door switches automatically shut off the UV light source should one of the doors be accidentally opened. An indicator light informs the user that cleaning is in process.



Part No.	Description	Unit
EIS-GGSN10	Goggle Sanitizer Cabinet	1

Kimberly-Clark Professional

Kimwipe® EX-L Delicate Task Wipers

- Soft & gentle on surfaces
- Low-linting



Part No.	Description	Pk	Cs
KCC-34155	Kimwipes Delicate Task Wipers, White, 4.4" x 8.4"	280	16,800
KCC-34133	Kimwipes Delicate Task Wipers, White, 11.8" x 11.8"	198	2,970
KCC-34256	Kimwipes Delicate Task Wipers, White, 14.7" x 16.6"	144	2,160

Haier

Pharmacy Refrigerators

- Each unit uses eco-friendly refrigerants and offer low power consumption relative to performance
- Features polyurethane insulation
- Comes with multiple alarm functions, including audible/visual alerts for monitoring, and some include remote alarm contacts
- Uses hermetically sealed compressors and forced-air circulation for efficient cooling.
- Adjustable wire shelves are standard across all models.
- Integrated digital control panels with microprocessors are used in all models, allowing for precise temperature control and easy monitoring
- All models are designed for quiet operation (35 – 43 dB(A))



Part No.	Description	Unit
HAI-HYC-1378	Pharmacy Refrigerator +2 to +8C, 1378L (49cf), Double glass door, 115V/60Hz	1
HAI-HYC-461FD	Pharmacy Refrigerator +2 to +8C, Energy Star, 461L (16.3cf), solid door, 115V/60Hz, 5 shelves, uniformity ±1.5C EA	1
HAI-HYC-940	Pharmacy Refrigerator +2 to +8C, 890L, Double glass door, 115V/60Hz	1

Labconco

FreeZone Freeze Dryers

- Intuitive 3.75" (9.5 cm) display with touch-pad control
- 4 languages: French, German, Spanish, English
- Swivel 180™ Rotating Controller
- On-board connectivity: USB-A & B, Ethernet, RS232/addressable RS485, and external temperature probe
- On-screen prompts
- Automatic performance optimization and specific heat tuning



Part No.	Description	Unit
LCC-710201000	FreeZone 2.5L -84° C Benchtop Freeze Dryer with stainless steel coil and collector 115V 60Hz	1
LCC-720401000	FreeZone 4.5L -105° C Benchtop Freeze Dryer with stainless steel coil and collector 115V 60Hz	1
LCC-780701000	FreeZone Small Tray Dryer 115V 50/60Hz	1
LCC-7368700	FreeZone Cart	1

United Scientific Supplies

Polypropylene Laboratory Tray

Multipurpose Polypropylene molded trays can be used for sterilizing, drying glassware, porcelain ware, and other laboratory supplies. These steam autoclavable trays have tapered walls designed for nesting.



Part No.	Description	Unit
USS-81701	Polypropylene Laboratory Tray, 18" x 14" x 3"	10
USS-81702	Polypropylene Laboratory Tray, 15" x 12" x 3"	10
USS-81731	Polypropylene Carrier Tray, 15" x 9.5" x 4.5"	6
USS-81721	Polypropylene Laboratory Tray, 15" x 14" x 5"	6
USS-LTRAY1	Polypropylene Laboratory Tray, Flat	10
USS-78103	Polypropylene Draining Tray, 16" x 16" x 4"	6

Hanna Instruments

Portable Meter Kit for Soil

Designed to bring simplicity to pH testing of soil, the HI99121 Direct Soil pH Portable Meter is designed specifically for soil analysis. The specialized probe features a rugged conical tip, open junction design, and a glass body making it ideal for measuring the pH of soil directly or after preparation of a soil slurry with deionized water.



Part No.	Description	Unit
HAN-HI99121	Direct Soil Measurement Ph Portable Meter	1

3M

Petrifilm™ Rapid Yeast and Mold Count Plate

- Utilizes convenience and efficiency of 3M™ Petrifilm™ plate platform
- Proven cost effective tool for detecting yeast and mold
- Can be used with both low and high water activity foods
- Eliminates need for labor intensive, multiple plate spread plate technique
- Versatility for air, swab or surface contact environmental sampling



Part No.	Description	Unit
NSS-FM-734-10	Yeast and Mold for 3M Petrifilm	10

J.G. Finneran

QuEChERS Food Extraction & Dispersive SPE Kits

- QuEChERS stands for (Quick, Easy, Cheap, Effective, Rugged & Safe)
- QuEChERS make it simple to prepare your food samples for analysis
- Extraction kits contain preweighed salts, so you can add them after the Acetonitrile Extraction step
- Dispersive SPE kits are assembled in 2ml and 15ml sizes



Part No.	Description	Unit
FIN-JG-F00115P	15 ml QuEChERS 1200mg MgSO ₄ , 400mg PSA Phase 2 Extraction Kit for AOAC 2007.01	50
FIN-JG-F00215P	15 ml QuEChERS 1200mg MgSO ₄ , 400mg PSA, 400mg C18 Phase 2 Extraction Kit for AOAC 2007.01	50
FIN-JG-F00315P	15 ml QuEChERS 1200mg MgSO ₄ , 400mg PSA, 400mg GCB Phase 2 Extraction Kit for AOAC 2007.01	50
FIN-JG-F00415P	15 ml QuEChERS 1200mg MgSO ₄ , 400mg PSA, 400mg C18, 400mg GCB Phase 2 Extraction Kit for AOAC 2007.01	50
FIN-JG-F00502P	2 ml QuEChERS 150mg MgSO ₄ , 50mg PSA Phase 2 Extraction Kit for AOAC 2007.01	100
FIN-JG-F00602P	2 ml QuEChERS 150mg MgSO ₄ , 50mg PSA, 50mg C18 Phase 2 Extraction Kit for AOAC 2007.01	100
FIN-JG-F00702P	2 ml QuEChERS 150mg MgSO ₄ , 50mg PSA, 50mg GCB Phase 2 Extraction Kit for	100
FIN-JG-F00715P	15mL QuEChERS Centrifuge Tube, 900mg MgSO ₄ , 150mg PSA, 45mg GCB	50
FIN-JG-F00802P	2 ml QuEChERS 150mg MgSO ₄ , 50mg PSA, 50mg C18, 50mg GCB Phase 2 Extraction Kit for AOAC 2007.01	100

Nalgene

PP Wide Mouth Laboratory Bottles

Nalgene™ Wide Mouth Leakproof Plastic Bottles are unique in design and very durable. They have a wide mouth making them very easy to fill. Bottles include a PP cap. Use only the cap provided to ensure leakproof performance.



Part No.	Description	Pk	Cs
NAL-2105-0004	Wide-Mouth Bottle PP 125 mL	12	72
NAL-2105-0008	Wide-Mouth Bottle PP 250 mL	12	72
NAL-2105-0032	Wide-Mouth Bottle PP 500 mL	12	72
NAL-2105-0032	Wide-Mouth Bottle PP 1000 mL	6	24
NAL-2121-0010	Wide-Mouth Bottle PP 4 L	-	6

GVS

Petri Dishes

- Non-cytotoxic virgin polystyrene, excellent optical clarity
- Engineered for optimum flatness to provide uniform agar thickness.
- Single or more vents are available: no-vent for anaerobic and long-term work, single vent for limiting gas exchange, multi vents for improving gas exchange.

Part No.	Description	Unit
GVS-PETRIDC340110380330	Petri Dishes Carrier Rack 340x110x380, 3 sections, hold 30 PD	5
GVS-PETRIDGC060151N-VS	Glass Petri Dishes, Cilindar, 60mm, 15mm, 1room, Sterile	288
GVS-PETRIDGC075151N-VS	Glass Petri Dishes, Cilindar, 75mm, 15mm, 1room, Sterile	-88
GVS-PETRIDGC090181N-VS	Glass Petri Dishes, Cilindar, 90mm, 18mm, 1room, Sterile	192
GVS-PETRIDGC100201N-VS	Glass Petri Dishes, Cilindar, 100mm, 20mm, 1room, Sterile	192
GVS-PETRIDPC035151N-VS	Petri Dishes, Cilindar, 35mm, 15mm, 1room, Sterile	3000
GVS-PETRIDPC049131N-VS	Petri Dishes, Cilindar, 49mm, 13mm, 1room, Sterile	1000
GVS-PETRIDPC055151N-VS	Petri Dishes, Cilindar, 55mm, 15mm, 1room, Sterile	1000
GVS-PETRIDP-C0601513VS	Petri Dishes, Cilindar, 60mm, 15mm, 1room, 3vents, Sterile	1000
GVS-PETRIDP-C0651513VS	Petri Dishes, Cilindar, 65mm, 15mm, 1room, 3vents, Sterile	1000
GVS-PETRIDPC070151N-VS	Petri Dishes, Cilindar, 70mm, 15mm, 1room, Sterile	800
GVS-PETRIDP-C0901513VS	Petri Dishes, Cilindar, 90mm, 15mm, 1room, 3vents, Sterile	500
GVS-PETRIDP-C0901523VS	Petri Dishes, Cilindar, 90mm, 15mm, 2room, 3vents, Sterile	500
GVS-PETRIDP-C0901533VS	Petri Dishes, Cilindar, 90mm, 15mm, 3room, 3vents, Sterile	500
GVS-PETRIDP-C1501513VS	Petri Dishes, Cilindar, 150mm, 15mm, 1room, 3vents, Sterile	200
GVS-PETRID-PCHW50090S	Petri Dishes, 90mm, 15mm, 1room, 3vents, Sterile	500
GVS-PETRIDPS130151N-VS	Petri Dishes, Square, SQ 130mm, 15mm, 1room, Sterile	150



Hanna Instruments

Humidity/Conductivity Meter

- The graphic LCD is easy to read outdoors in bright sunlight as well as in low-lit areas with the use of the backlight. A combination of dedicated and virtual keys allows easy, intuitive operation in a choice of languages.
- The compact, durable carrying case is thermoformed to hold all necessary components for taking a field measurement, including the meter and electrode, beakers, buffer solutions, cleaning solutions, and storage solutions.



Part No.	Description	Unit
HAN-HI98190	Professional Waterproof Portable pH/ORP Meter	1
HAN-HI98191	Professional Waterproof Portable pH/ORP/ISE Meter	1
HAN-HI98192	Waterproof portable EC/resistivity/TDS/NaCl auto ranging meter with graphic display; range	1

Rotronic

Water Activity Meter

A hand-held water activity indicator with data logging capability. The HP23-AW-A is useful in a number of applications such as spot check measurements to confirm product stability and indicate shelf life.



Part No.	Description	Unit
JAD-HP23-AW-A	Water Activity Meter	1

GVS

Bottle-Top Filters- ZapCap™

GVS Bottle-top Filters is ideal solution for the filtration of cell culture media and HPLC media solution. ZapCap™ filters are equipped with side tubing nozzle (bottle-top). This ready to use filter is available with prefilter too.



Part No.	Description	Unit
GVS-10443401	ZapCap-S, Cellulose Acetate, 0.2µm, STERILE, Polystyrene	12
GVS-10443411	ZapCap-S, Cellulose Acetate, 0.45µm, STERILE, Polystyrene-S	12
GVS-10443421	ZapCap-CR, PA6.6, 0.2µm, NON STERILE, Polypropylene	12
GVS-10443423	ZapCap-CR, PA6.6, 0.45µm, NON STERILE, Polypropylene	12
GVS-10443425	ZapCap-CR, PTFE, 0.45µm, NON STERILE, Polypropylene	12
GVS-10443430	ZapCap-S PLUS, Cellulose Acetate, 0.2µm, STERILE, Cellulose Acetate with Glass Fiber Prefilter	12
GVS-10443435	ZapCap-S PLUS, Cellulose Acetate, 0.45µm, STERILE, Cellulose Acetate with Glass Fiber Prefilter	12

BrandTech® Dispensette® S Digital Bottletop Dispensers

The Dispensette® S builds on a half century of BRAND® bottletop dispenser experience. It features a floating piston design with even lower operating forces, a new safety discharge system, and an improved sealing system, eliminating sealing washers for even better reliability. Designed to attach directly to your reagent bottle eliminating poured transfer. In addition to the product improvements, a new 0.1-1mL size is available.




Part No.	Description	Unit
BRD-4600310	Dispensette S, Digital w/standard valve, 0.1-1ml	1
BRD-4600320	Dispensette S, Digital w/standard valve, 0.2-2ml	1
BRD-4600330	Dispensette S, Digital w/standard valve, 0.5-5ml	1
BRD-4600340	Dispensette S, Digital w/standard valve, 1-10ml	1
BRD-4600350	Dispensette S, Digital w/standard valve, 5-50ml	1
BRD-4600311	Dispensette S, Digital w/recirculation valve, 0.1-1ml	1
BRD-4600321	Dispensette S, Digital w/recirculation valve, 0.2-2ml	1
BRD-4600331	Dispensette S, Digital w/recirculation valve, 0.5-5ml	1
BRD-4600341	Dispensette S, Digital w/recirculation valve, 1-10ml	1
BRD-4600351	Dispensette S, Digital w/recirculation valve, 2.5-25ml	1
BRD-4600361	Dispensette S, Digital w/recirculation valve, 5-50ml	1

Sol-Millennium Syringe without Needle

- High-transparency barrel to easily visualize contents
- Silicone lubrication ensures smooth plunger movement
- Plunger backstop that prevents accidental withdrawal of plunger from the barrel
- Available in Luer Lock, Luer Slip Tip, and Catheter Tip options

Part No.	Description	Unit
SOL-P180001	1ml Luer Lock Syringe without needle	800
SOL-P180003	3ml Luer Lock Syringe without needle	800
SOL-P180005	5ml Luer Lock Syringe without needle	800
SOL-P180010	10ml Luer Lock Syringe without needle	1,200
SOL-P180020	20ml Luer Lock Syringe without needle	800
SOL-P180030	30ml Luer Lock Syringe without needle	800
SOL-P180060	60ml Luer Lock Syringe without needle	800



United Scientific Supplies Magnetic Bar Retriever

350mm length magnetic bar retriever.

Part No.	Description	Unit
USS-MBR350	Magnetic Bar Retriever, 350mm	1



Globe Scientific Supplies Stir Bar Retriever

Easily recover magnetic stir bars from sensitive or corrosive solutions with Globe Scientific's Stir Bar Retriever. This essential lab tool features a powerful magnet permanently embedded in a chemically resistant PTFE (polytetrafluoroethylene) rod, ensuring durability, safety, and compatibility with a wide range of laboratory liquids.



Part No.	Description	Unit
GLO-SWN652	Stir Bar, Cylindrical, PTFE, 6mm	2
GLO-SWN654	Stir Bar, Cylindrical, PTFE, 8mm	2
GLO-SWN656	Stir Bar, Cylindrical, PTFE, 12mm	2
GLO-SWN658	Stir Bar, Cylindrical, PTFE, 15mm	2
GLO-SWN660	Stir Bar, Cylindrical, PTFE, 20mm	2
GLO-SWN662	Stir Bar, Cylindrical, PTFE, 25mm	2
GLO-SWN664	Stir Bar, Cylindrical, PTFE, 30mm	2
GLO-SWN666	Stir Bar, Cylindrical, PTFE, 35mm	2
GLO-SWN668	Stir Bar, Cylindrical, PTFE, 40mm	2
GLO-SWN670	Stir Bar, Cylindrical, PTFE, 45mm	2
GLO-SWN672	Stir Bar, Cylindrical, PTFE, 50mm	2
GLO-SWN674	Stir Bar, Cylindrical, PTFE, 60mm	2
GLO-SWN676	Stir Bar, Cylindrical, PTFE, 70mm	2

GVS Syringe Filter PTFE

Part No.	Description	Unit
GVS-FJ25BNP-PH004AD01	ABLUO Syringe Filter NON STERILE 25mm PTFE Hydrophilic, 0.45µm	500
GVS-FJ25B-NPPT002AD01	ABLUO Syringe Filter NON STERILE 25mm PTFE Hydrophobic, 0.2µm	500
GVS-FJ25B-SCPV004AL01	ABLUO Syringe Filter STERILE 25mm PVDF, 0.45µm	1000
GVS-FJ25B-SCPV002AL01	ABLUO Syringe Filter STERILE 25mm PVDF, 0.2µm	50



Ahlstrom-Munksjo Qualitative Filter Papers, Ahlstrom 631

- Quantitative analysis
- Low ash content 0.06%
- Untreated
- High wet-strength



Part No.	Description	Unit
AHL-6150-0550	Qualitative Filter Paper -5.5cm	100
AHL-6150-0700	Qualitative Filter Paper -7cm	100
AHL-6150-0750	Qualitative Filter Paper -7.5cm	100
AHL-6150-0900	Qualitative Filter Paper -9cm	100
AHL-6150-1250	Qualitative Filter Paper -12.5cm	100
AHL-6150-1500	Qualitative Filter Paper -15cm	100
AHL-6150-1850	Qualitative Filter Paper -18.5cm	100
AHL-6150-2050	Qualitative Filter Paper -20.5cm	100
AHL-6150-2400	Qualitative Filter Paper -24cm	100
AHL-6150-2500	Qualitative Filter Paper -25cm	100

Brandtech Scientific 96-Well PCR Plates, Low-Profile

- Clear or White
- 96 x 0.1ml wells (low-profile)
- Flat top for ideal sealing with automated & manual seals
- Optically clear with black printed alphanumeric
- Certified DNase, RNase, and Pyrogen Free

Part No.	Description	Unit
BRD-781377	96-well PCR plate skirt Low Profile, Clear	500
BRD-781378	96-well PCR plate skirt Low Profile, White	500
BRD-781378	96-well PCR plate skirted Low Profile	500

MP Biomedicals FastPrep-24™ Classic

- Simultaneous homogenization of up to 24 samples within 40 seconds
- Thorough and quick lysis of any tissues and cells
- 9 interchangeable adapters allows a unique versatility in sample number, size and temperature conditions



Part No.	Description	Unit
MPB-116004500	FastPrep-24 Classic bead beating grinder and lysis system	1

American Scientific Products Low Retention Pipette Tips

- Non-Sterile options available
- Universal fit on research grade pipettors
- Accurate, precise and reliable
- Certified free of RNase, DNase, DNA, and pyrogens
- Extended length tips enable users to reach the bottom of tubes without contaminating pipettors
- Aerosol barrier will not trap liquid or contaminate samples in the event of over-pipetting.
- Graduated tips provide visible reference points

Part No.	Description	Unit
ASP-MT-10XL-96	MoleBio 10ul LongReach Universal Filter Pipet Tips, PP, Sterile	960
ASP-MT-20-96	MoleBio 20ul Universal Filter Pipet Tips, PP, Sterile	960
ASP-MT-100-96	MoleBio 100ul Universal Filter Pipet Tips, PP, Sterile	960
ASP-MT-200-96	MoleBio 200ul Universal Filter Pipet Tips, PP, Sterile	960
ASP-MT-300-96	MoleBio 300ul Universal Filter Pipet Tips, PP, Sterile	960
ASP-MT-1250-96	MoleBio 1250ul Universal Filter Pipet Tips, PP, Sterile	960



American Scientific Products 15ml & 50ml Polypropylene Centrifuge Tubes

Centrifuge Tubes are ideal for transporting and collecting specimens. Manufactured from chemically inert materials, each tube is suitable as the primary containers for shipping biological specimens. The screw caps feature a continuous 1½-turn thread that assist in securing the cap while keeping it easy to open and close. A molded-in sealing ring wraps the tube wall to prevent liquid bypass and provides a positive stop for the cap. Not autoclavable.

Part No.	Description	Unit
ASP-3216P	Centrifuge Tubes, 15 ml, Screw Cap	500
ASP-3217	Centrifuge Tube, 15 ml, Molded in Graduations, Sterile, Case	500
ASP-50CTMB	50ml Centrifuge Tube, Non-Sterile, Molded Grads., w/ screw cap	500
ASP-50CTMBS	50ml Centrifuge Tube, Sterile, Molded Grads., w/ screw cap	500
ASP-50CTPB	50ml Centrifuge Tube, Non-Sterile, Printed White Grads., w/ screw cap	500
ASP-50CTPBS	50ml Centrifuge Tube, Sterile, Printed White Grads., w/ screw cap	500
ASP-50CTSSMB	50ml Centrifuge Tube, Self Standing, Non-Sterile, Molded Grads., w/ screw cap	500



Corning Volumetric Flask

Part No.	Description	Unit
COR-5580-2	PYREX 2 mL Class A Volumetric Flask	12
COR-5640-10FO	PYREX 10 mL Class A Volumetric Flask	12
COR-5640-50FO	PYREX 50 mL Class A Volumetric Flask	1
COR-5640-100FO	PYREX 100 mL Class A Volumetric Flask	1
COR-5640-250FO	PYREX 250 mL Class A Volumetric Flask	1
COR-5640-1LFO	PYREX 1L Class A Volumetric Flask	1



J.G. Finneran Wide Opening Glass Vials

This 11mm crimp vial is made of amber borosilicate glass with White Graduated Spot, with a 40% larger opening that prevents broken needles due to increased target area. Designed to work in Agilent and other autosamplers. Sample capacity of 2mL, 12mm in diameter, and 32mm in height. Vials compatible with 11mm crimp and Poly Crimp™ Seals. 1,000 vials per case with 10 shelf packs of 100 vials packed in plastic shrink wrapped scallop packaging.



Part No.	Description	Unit
FIN-3211LE-1232A	2.0mL LO (Large Opening) Amber Vial, 12x32mm, with White Graduated Spot, 11mm	1000
FIN-3211LE-12AZ	Silanized - 2.0mL LO (Large Opening) Amber Vial, 12x32m	1000

J.G. Finneran 20mL Headspace Vials

This 11mm crimp vial is made of amber borosilicate glass with White Graduated Spot, with a 40% larger opening that prevents broken needles due to increased target area. Designed to work in Agilent and other autosamplers. Sample capacity of 2mL, 12mm in diameter, and 32mm in height. Vials compatible with 11mm crimp and Poly Crimp™ Seals. 1,000 vials per case with 10 shelf packs of 100 vials packed in plastic shrink wrapped scallop packaging.



Part No.	Description	Unit
FIN-320020-2375	20mL Clear Headspace Vial, 23x75mm, Flat Bottom, 20mm Beveled Crimp Top	1000
FIN-322AXLN-2375	20mL Clear Headspace Vial, 23x75mm (for CTC PAL, Perkin Elmer), Beveled Bottom, 20mm Flat Top Crimp, Long Neck	1000

J.G. Finneran Pre-Assembled R.A.M.™ Vials with 9mm Caps

Components packaged in a tray that fits into bench drawer for easy access. Packaging keeps products visible & particle free. Pre-assembled. Packs come ready-to-use.



Part No.	Description	Unit
FIN-80494-12F	Clear 9mm Thread R.A.M.™ Vial, 12x32mm, 0.040" PTFE/Butyl Rubber Septa, Ribbed Black Cap	100
FIN-80494E-12F	Clear 9mm Thread R.A.M.™ Vial w/Grad. Spot, 12x32mm, PTFE/Butyl Rubber Septa, Ribbed Black Cap	12

DWK Life Sciences Distilling Heads

- With a Condenser and Medium Length Standard Taper
- Joints Standard Taper 14/20 outer joint on the vertical tube can be converted to a thermometer joint using the 273500, size 21 bushing adapter.
- The height of this joint is fixed in relation to the side arm, permitting use of 25 mm immersion thermometers 14/35 Condenser Joint.
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements.



Part No.	Description	Unit
KMC-286900-0000	Kimble 286900-0000 Kontes Distillation Head with Condenser and Medium Length 14/20 [ST] Joints	1

Dynalab Corp Drying Rack

- Small (3"), medium (4") and large (5") pegs can be moved and repositioned to hold a variety of labware including beakers, flasks, bottles, test tubes and cylinders
- Rugged ABS plastic will never corrode and easily wipes clean
- Drip/storage trays have scalloped wells to hold long thin labware
- Extra accessories and a variety of different sizes/configurations let you create a drying system that works for you



Part No.	Description	Unit
DYN-259184	Drying Rack, 72 place	1

United Scientific Supplies Low Form Beakers

Heavy duty beakers are ideal for tough handling conditions including mechanized washing. Beakers feature a double graduated metric scale and a marking area.



Part No.	Description	Unit
USS-BG1000-100	Beaker, 100 ml, Low Form, Borosilicate Glass	12
USS-BG1000-250	Beaker, 250 ml, Low Form, Borosilicate Glass	12
USS-BG1000-500	Beaker, 500 ml, Low Form, Borosilicate Glass	6
USS-BG1000-3000	Beaker, 3000 ml, Low Form, Borosilicate Glass	1

Globe Scientific Wide Mouth Sample Bottles

These rupture-resistant bottles are extremely rigid and feature a wide-mouth opening for easy filling and emptying of liquids and powders. The deep threaded neck allows the screw cap to seal on tight, offering exceptional leak-resistance especially during shipping. This bottle is approved for use with food products.



Part No.	Description	Unit
GLO-601623	Bottle w/ Screwcap, Wide Mouth, PP, Graduated, 100mL	1
GLO-601625	Bottle w/ Screwcap, Wide Mouth, PP, Graduated, 250mL	1
GLO-601627	Bottle w/ Screwcap, Wide Mouth, PP, Graduated, 500mL	1
GLO-601631	Bottle w/ Screwcap, Wide Mouth, PP, Graduated, 2000mL	1

CELLTREAT Single Channel Pipette

- High precision adjustable pipette with easy-to-read volume display
- Continuously adjustable with locking ring to prevent volume changes while pipetting
- Lightweight design and easy operation with either hand
- Reduced operating forces; plunger is easy to press and dispense
- Stainless steel, PTFE, and PVDF components for long life
- Autoclavable up to 6 times



Part No.	Description	Unit
CEL-230299	1-10µL Pipette, Adjustable Volume, Single Channel	1
CEL-230305	20-200µL Pipette, Adjustable Volume, Single Channel	1
CEL-230310	100-1000µL Pipette, Adjustable Volume, Single Channel	1

United Scientific Supplies Pasteur Pipettes, LDPE

Heavy duty beakers are ideal for tough handling conditions including mechanized washing. Beakers feature a double graduated metric scale and a marking area.



Part No.	Description	Unit
USS-P31201	Pasteur Pipettes, LDPE, 1.0mL, 150mm	500
USS-P31205	Pasteur Pipettes, LDPE, 3.0mL, 150mm	500

Bel-Art Products, Inc. Plastic Büchner Funnels

Heavy duty beakers are ideal for tough handling conditions including mechanized washing. Beakers feature a double graduated metric scale and a marking area.



Part No.	Description	Unit
BEL-14676-0043	Funnel, PP, Standard Stem, 40mm	6
BEL-14676-0100	Funnel, PP, Standard Stem, 100mm	4
BEL-14608-0000	Funnel, HDPE, Standard Stem, 50mm	1
BEL-14609-0000	Funnel, HDPE, Standard Stem, 150mm	1
BEL-14911-0000	Funnel, HDPE, Standard Stem, 1000mm	1

Corning PYREX® Brand 5340 Filtering Flasks

These PYREX® flasks have sidearm tubulations and are blown in special molds designed to give maximum mechanical strength. For convenience, the flasks are graduated to show approximate capacity. The neck finish affords an excellent fit for rubber stoppers.



Part No.	Description	Unit
BEL-14676-0043	Funnel, PP, Standard Stem, 40mm	6
BEL-14676-0100	Funnel, PP, Standard Stem, 100mm	4
BEL-14608-0000	Funnel, HDPE, Standard Stem, 50mm	1
BEL-14609-0000	Funnel, HDPE, Standard Stem, 150mm	1
BEL-14911-0000	Funnel, HDPE, Standard Stem, 1000mm	1

SPEX CertiPrep

21 Component Terpene Mix in Methanol

SPEX CertiPrep offers ISO/IEC 17025 and ISO Guide 34 Certified Reference Materials (CRMs) for all of the common contaminants such as pesticide residues, residual solvents, and heavy metals, as well as qualitative analysis CRMs, such as terpenes.

Part No.	Description	Size	Unit
SPX-CAN-TERP-MIX1H	TERPENES MIX 1- HIGH Terpenes Mix 1 High Level 21 Component Mix	1000ug/ml	1
SPX-CAN-TERP-MIX2H	TERPENES MIX 1- HIGH Terpenes Mix 1 High Level 21 Component Mix	1000ug/ml	1

Ethyl Alcohol, 200 Proof

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-D9080TP	Ethyl Alcohol 200PF, Poly Bottle	1 Gal	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-EX0276-3	Ethyl Alcohol, 200 Proof GR ACS	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-AH090-4	Ethyl Alcohol, 200 Proof	4L	1

n-Pentane, 99%

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-E9880	Pentane, N-, Certified	5 Gal	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-34956-4L	Pentane suitable for HPLC, ≥99.0%	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-34956-1L	Pentane CHROMASOLV®, for HPLC, >99.0%	1L	1

Acetonitrile

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-C0157	Acetonitrile, 99.9% Min, HPLC/ACS	4L	1
CP-C9725-EA	Acetonitrile, LC-MS & HPLC, >99.9%	4L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-AX0145-1	Acetonitrile HPLC Grade	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-34851-4L	Acetonitrile CHROMASOLV, for HPLC, gradient grade, >99.9%	4L	1

PanReac

Part No.	Description	Size	Unit
PRA-701881.1611	Acetonitrile for LC-MS	1L	1
PRA-701881.16153	Acetonitrile for LC-MS	4L	4
PRA-701881.1612	Acetonitrile for LC-MS	2.5	1

Methanol

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-C2305EA	Methanol, >99.9%, HPLC	4L	1
CP-C7142	Methanol, 99.8% Min, ACS	4L	1
CP-E7142	Methanol, ACS, Metal Pail	19L	1
CP-T7142	Methanol, ACS	55 Gal	1
CP-C9635	Methanol, UHPLC/LC-MS	4L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-MX0486-1	Methanol, OmniSolv(R) LC-MS	4L	1
SIG-34885-1L-R	Methanol, HPLC, Gradient Grade, ≥99.9%	1L	1

Honeywell

Part No.	Description	Size	Unit
SIG-MX0486-1	Methanol, LC-MS	4L	1
HBJ-34966-4L	Methanol, LC-MS Chromasolv(R)	4L	1
HBJ-230-4-CUS04878	Methanol, HPLC	4L	1

PanReac

Part No.	Description	Size	Unit
PRA-701091.1611	Methanol for LC-MS	1L	1
PRA-701091.1612	Methanol for LC-MS	4L	4
PRA-701091.1613	Methanol for LC-MS	2.5	1

Isopropyl Alcohol

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-T9091	Isopropyl Alcohol >99.5% Meets USP, Laboratory Reagent	1L	1
CP-E9093	Isopropyl Alcohol, ACS, >99%, Metal Pail	4L	1
CP-C9093	Isopropyl Alcohol, ACS, >99%	4L	1
CP-A9652-CS	Isopropyl Alcohol, ACS/USP/NF, >99%	1 Pint	24

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-P9416-100ML	Isopropyl Alcohol, 70% Reagent	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-323-1L	Isopropyl Alcohol, HPLC	1L	1
HBJ-323-4-CUS04878	Isopropyl Alcohol, High Purity 99.9%	4L	4

2-Propanol, 99.9%

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-34863	2-Propanol suitable for HPLC, 99.9%	4L	1
SIG-650447	2-Propanol CHROMASOLV® Plus, for HPLC, 99.9%	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-34863-4L	2-Propanol CHROMASOLV®, 99.9%	4L	1
HBJ-34965-4L	2-Propanol LC-MS CHROMASOLV®	4L	1
HBJ-650447-4L	2-Propanol CHROMASOLV® Plus, for HPLC, 99.9%	4L	1

Acetone

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-C9043	Acetone, 99.9% Min, HPLC/ACS	4L	1
CP-E7042	Acetone, ACS, >99.5%	19L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-010-4	Acetone for HPLC, GC	4L	1
HBJ-650501-1L	Acetone CHROMASOLV® Plus, for HPLC, >99.9%	1L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-179124-500ML	Acetone ACS reagent, ≥99.5%	500ML	1
SIG-179973-4L	Acetone Laboratory Reagent, ≥99.5%	4L	1

Sodium Hydroxide

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-B2435P-EA	Sodium Hydroxide 0.1N (0.1M), 1L	1L	1
CP-C2435P-EA	Sodium Hydroxide, 0.1N 4L	4L	1
CP-E2435-EA	Sodium Hydroxide 0.1N (0.1M), 20 L	20L	1
CP-A2438P	Sodium Hydroxide 1.0N (1.0M) 500ml	500mL	1
CP-C2438P-EA	Sodium Hydroxide 1.0N (1.0M) 4L	4L	1
CP-E2438	Sodium Hydroxide 1.0N (1.0M) 20L	20L	1
CP-E2393-EA	Sodium Hydroxide 10N, 20L	20L	1
CP-A2390D	Sodium Hydroxide, Pellets, >97.0%, ACS, 500g	500g	1
CP-M2390D	Sodium Hydroxide, Pellets, >97.0%, ACS, 2.5kg	2.5Kg	1
CP-I2390D	Sodium Hydroxide, Pellets, >97.0%, ACS, 12kg	12Kg	1

Ethanol 200 Proof CDA Formula 12A

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-T9901	Ethanol CDA12A 200 Proof (N-Heptane) High Purity Extraction Grade	55 Gal	1

Ethanol (190 & 200 Proof)

Volu-Sol

Part No.	Description	Size	Unit
VOL-V95-640	Ethanol, 190 Proof 95% non-denatured	5 Gal	1
VOL-V4B-005	Pure Denatured SDA Ethanol (200 Proof), Formula 40-B	55 Gal	1
VOL-VET-055TK	Pure Ethanol Absolute, 200 Proof, Tax-Paid ACS, USP, FCC, Koshe	55 Gal	1

**Cases Available, Bulk
Quantity & Discounts
Provided**

TWEEN 20**ChemPure Chemicals**

Part No.	Description	Size	Unit
CPB-P20370-1L	Tween 20	1L	1
CPB-P20370-4L	Tween 20	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-66117	Tween 20	4L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-P9416-100ML	Tween 20	100ML	1
SIG-P9416-250ML	Tween 20	250ML	1

Formic Acid**ChemPure Chemicals**

Part No.	Description	Size	Unit
CP-V9623	Formic Acid, LC-MS Grade, 10x1mL Ampoule	1mL	10
CP-L9636	Formic Acid, LC-MS	50mL	1
CP-B0615	Formic Acid 98-100% High Purity, ACS	1L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-695076-100ML	Formic acid, ACS reagent, ≥96%	100mL	1
SIG-695076-500ML	Formic acid, ACS reagent, ≥96%	500mL	1

Honeywell

Part No.	Description	Size	Unit
HBJ-56302-50ML-GL	Formic Acid, LC-MS	50mL	1
HBJ-56302-1L-GL	Formic Acid, LC-MS	1L	1

Nitric Acid**ChemPure Chemicals**

Part No.	Description	Size	Unit
CP-M1778P	Nitric Acid, 70% ACS/High Purity/Trace, Polybottle	2.5L	1
CP-M1780P	Nitric Acid, 70% ChemTrace, Polybottle	2.5L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-438073	Nitric Acid ACS reagent, 70%	2.5L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-438073-2.2L-P	Nitric Acid ACS reagent, 70%	2.2L	1
HBJ-84385-2.5L	Nitric acid for Trace Analysis	2.5L	1

Hydrochloric Acid**Sigma-Aldrich**

Part No.	Description	Size	Unit
SIG-258148-2.5L-GL	Hydrochloric acid ACS reagent, 37%	2.5L	1
SIG-258148-4L	Hydrochloric acid ACS reagent, 37%	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-08256-1L	Hydrochloric acid >30%, TraceSELECT®, for trace analysis	1L	1
HBJ-08256-2.5L	Hydrochloric acid >30%, TraceSELECT®, for trace analysis	2.5L	1

ChemPure Chemicals

Part No.	Description	Size	Unit
CP-D1494P-EA	Hydrochloric Acid, ACS/High Purity/Trace	3.84L	1
CP-M1498P	Hydrochloric Acid, ChemTrace grade, impurities at the ppb level	2.5L	1

Water**ChemPure Chemicals**

Part No.	Description	Size	Unit
CP-C0001	Water, HPLC	4L	1
CP-C9605-EA	Water, UHPLC/LC-MS for Mass Spectrometry	4L	1

Sigma-Aldrich

Part No.	Description	Size	Unit
SIG-270733	Water, HPLC	4L	1
SIG-900682	Water, UHPLC/LC-MS for Mass Spectrometry	4L	1

Honeywell

Part No.	Description	Size	Unit
HBJ-365-4-EA	Water, HPLC	4L	1
HBJ-365-4-CUS04878	Water, UHPLC/LC-MS for Mass Spectrometry	4L	1

PanReac

Part No.	Description	Size	Unit
PRA-701074.1611	Water for LC-MS	1L	1
PRA-701074.1612	Water for LC-MS	2.5L	1
PRA-701074.1613	Water for LC-MS	4L	4

PanReac
AppliChem
ITW Reagents

Jade Scientific is the exclusive distributor of PanReac LC-MS Grade Chemicals in the United States.

Laboratory Furniture

Jade Scientific, Inc. offers a comprehensive selection of high-quality laboratory furniture tailored to the needs of cannabis manufacturing and testing facilities. Our durable, chemical-resistant benches, cabinetry, and workstations are designed to support safe, efficient workflows in both production and analytical environments. From modular casework to custom lab layouts, we provide flexible solutions that meet industry standards and adapt to evolving compliance requirements.

- Flexible & Modular Lab Furniture Systems
- Pre-plumbed & Wired Workstations
- Chemical Fume Hoods
- Shelving Systems
- Tables/Tabling Systems
- Healthcare
- Stainless Steel Product Line
- Work Surfaces



Cannabinoid Potency Standards

Jade Scientific, Inc. offers a comprehensive selection of Cannabinoid Potency Standards to support accurate and reliable cannabis testing. These high-purity reference materials are essential for quantifying key cannabinoids like THC, CBD, CBG, and more. Designed for use in analytical methods such as HPLC and GC, our standards ensure compliance, consistency, and confidence in every result.



Single Component CRMs:

Cannabinoid	CAS Number	Cayman Chemicals	Unit	Cerilliant	Unit	LGC	Unit
CBD	13956-29-1	CAY-ISO60156-1MG	1,000 µg/ml	C-045-1ML	1.0 mg/mL	DRE-A10946000ME-1000	1,000 µg/ml
CBD-A	1244-58-2	CAY-18090-1MG	1,000 µg/ml	C-144-1ML	1.0 mg/mL	DRE-A10946020AL-1000	1,000 µg/ml
CBD-V	24274-48-4	CAY-20165-1MG	1,000 µg/ml	C-140-1ML	1.0 mg/mL	DRE-10946040ME-1000	1,000 µg/ml
CBDV-A	31932-13-5	CAY-30245-1MG	1,000 µg/ml	C-152-1ML	1.0 mg/mL	DRE-A10946035AL-1000	1,000 µg/ml
CBG	25654-31-3	CAY-20164-1MG	1,000 µg/ml	C-141-1ML	1.0 mg/mL	DRE-A10946100ME-1000	1,000 µg/ml
CBG-A	25555-57-1	CAY-20019-1MG	1,000 µg/ml	C-142-1ML	1.0 mg/mL	DRE-A10946120AL-1000	1,000 µg/ml
CBN	521-35-7	CAY-ISO60183-1MG	1,000 µg/ml	C-046-1ML	1.0 mg/mL	DRE-A10946200ME-1000	1,000 µg/ml
CBN-A	2808-39-1	CAY-21734-1MG	1,000 µg/ml	C-153-1ML	1.0 mg/mL	DRE-A10946220ME-1000	1,000 µg/ml
THC	1972-08-3	CAY-ISO60157-1MG	1,000 µg/ml	T-005-1ML	1.0 mg/mL	DRE-A17405100ME-1000	1,000 µg/ml
THC-A	23978-85-0	CAY-ISO60175-1MG	1,000 µg/ml	T-093-1ML	1.0 mg/mL	DRE-A17405150AL-1000	1,000 µg/ml
THC-V	31262-37-0	CAY-18091-1MG	1,000 µg/ml	T-094-1ML	1.0 mg/mL	DRE-A17405170ME-1000	1,000 µg/ml
THCVA	39986-26-0	CAY-21259-1MG	1,000 µg/ml	T-111-1ML	1.0 mg/mL	DRE-A17405190AL-1000	1,000 µg/ml
D8-THC	5957-75-5	CAY-ISO60158-1MG	1,000 µg/ml	T-032-1ML	1.0 mg/mL	DRE-A17405050ME-1000	1,000 µg/ml
CBC	20675-51-8	CAY-ISO60163-1MG	1,000 µg/ml	C-143-1ML	1.0 mg/mL	DRE-A10945900ME-1000	1,000 µg/ml
CBC-A	185505-15-1	CAY-37558-1MG	1,000 µg/ml	C-244-1ML	1.0 mg/mL	DRE-A10945910AL-1000	1,000 µg/ml
CBC-V	24274-48-4	CAY-21974-1MG	1,000 µg/ml	NA	NA	DRE-A10946040ME-1000	1,000 µg/ml
CBCV-A	1628112-69-5	CAY-32718-1MG	1,000 µg/ml	C-256-1ML	1.0 mg/mL	DRE-A10945927AL-1000	1,000 µg/ml
CBL	21366-63-2	CAY-22036-1MG	1,000 µg/ml	C-154-1ML	1.0 mg/mL	DRE-A10945950AL-1000	1,000 µg/ml
CBC-O	55824-09-4	CAY-21742-1MG	1,000 µg/ml	C-251-1ML	1.0 mg/mL	NA	NA
CBL-A	40524-99-0	CAY-35745-1MG	1,000 µg/ml	C-264-1ML	0.5 µg/mL	DRE-A10945960AL-1000	1,000 µg/ml
CBT	31508-71-1	CAY-21295-1MG	1,000 µg/ml	C-241-1ML	1.0 mg/mL	DRE-A10945930AL-1000	1,000 µg/ml
(-)-Δ9-THC	1972-08-3	CAY-ISO60157-1MG	1,000 µg/ml	T-005-1ML	1.0 mg/mL	DRE-A17405100ME-1000	1,000 µg/ml
(-)-Δ9-THC-d3	81586-39-2	CAY-19332-1MG	1,000 µg/ml	T-003-1ML	100 µg/ml	NA	NA
CBD-d9	1246819-21-5	CAY-16203-1MG	1,000 µg/ml	NA	NA	DRE-A10946005AL-100	100 µg/ml
CBC-d9	Not Available	CAY-21294-100UG	100 µg/ml	NA	NA	NA	NA
CBD-d3	Not Available	CAY-10009523-100UG	100 µg/ml	C-084-1ML	100 µg/ml	NA	NA
CBN-d3	1435934-54-5	NA	NA	C-115-1ML	100 µg/ml	NA	NA

Cannabinoid Potency Standards

Single Component CRMs:

Cannabinoid	CAS Number	Restek	Unit	Spex/NSI	Unit	Absolute Standards	Unit
CBD	13956-29-1	RES-34011	1,000 µg/ml	S-10241	1,000 µg/ml	AES-73467	1,000 µg/ml
CBD-A	1244-58-2	RES-34094	1,000 µg/ml	S-10249	1,000 µg/ml	AES-73683	1,000 µg/ml
CBD-V	24274-48-4	RES-34123	1,000 µg/ml	S-10245	1,000 µg/ml	AES-74129	1,000 µg/ml
CBDV-A	31932-13-5	RES-34134	1,000 µg/ml	C-2516	100 µg/ml	AES-65267	100 µg/ml
CBG	25654-31-3	RES-34091	1,000 µg/ml	S-10246	1,000 µg/ml	AES-74038	1,000 µg/ml
CBG-A	25555-57-1	RES-34135	1,000 µg/ml	S-10247	1,000 µg/ml	AES-74130	1,000 µg/ml
CBN	521-35-7	RES-34010	1,000 µg/ml	S-10242	1,000 µg/ml	AES-73769	1,000 µg/ml
CBN-A	2808-39-1	RES-34138	1,000 µg/ml	C-2518	100 µg/ml	NA	NA
THC	1972-08-3	RES-34067	1,000 µg/ml	S-10260	1,000 µg/ml	AES-73040	1,000 µg/ml
THC-A	23978-85-0	RES-34111	1,000 µg/ml	S-11056	1,000 µg/ml	AES-73765	1,000 µg/ml
THC-V	31262-37-0	RES-34100	1,000 µg/ml	S-11057	1,000 µg/ml	AES-64811	1,000 µg/ml
THCVA	39986-26-0	RES-34139	1,000 µg/ml	C-2519	100 µg/ml	NA	NA
D8-THC	5957-75-5	RES-34090	1,000 µg/ml	S-10261	1,000 µg/ml	AES-73039	1,000 µg/ml
CBC	20675-51-8	RES-34092	1,000 µg/ml	S-10248	1,000 µg/ml	AES-73604	1,000 µg/ml
CBC-A	185505-15-1	RES-34137	1,000 µg/ml	NA	NA	AES-67491	100 µg/ml
CBC-V	24274-48-4	NA	NA	NA	NA	AES-74129	1,000 µg/ml
CBCV-A	1628112-69-5	NA	NA	NA	NA	NA	NA
CBL	21366-63-2	RES-34130	1,000 µg/ml	C-2517	100 µg/ml	NA	NA
CBC-O	55824-09-4	NA	NA	NA	NA	NA	NA
CBL-A	40524-99-0	RES-34136	500 µg/ml	NA	NA	NA	NA
CBT	31508-71-1	NA	NA	NA	NA	NA	NA
(-)-Δ9-THC	1972-08-3	RES-34067	1,000 µg/ml	S-10260	1,000 µg/ml	AES-73040	1,000 µg/ml
(-)-Δ9-THC-d3	81586-39-2	NA	NA	NA	NA	AES-96479	100 µg/ml
CBD-d9	1246819-21-5	NA	NA	NA	NA	NA	NA
CBC-d9	Not Available	NA	NA	NA	NA	NA	NA
CBD-d3	Not Available	NA	NA	NA	NA	NA	NA
CBN-d3	1435934-54-5	NA	NA	NA	NA	NA	NA



Cannabinoid Potency Standards

Multi Component CRMs:

		2 Part Mix				3 Part Mix			
		Absolute Standards	Cerilliant	LGC	LGC	Absolute Standards	Cayman	Cerilliant	
Cannabinoid	CAS Number	AES-99361A	CER-C-217-1ML	LGC-DRE-A50000199AL	LGC-DRE-A50000200AL	AES-97837	CAY-23251-1EA	CER-T-108-0.5ML	CER-C-220-1ML
CBD	13956-29-1		1.0 mg/mL	1,000 µg/ml		100 µg/ml	1,000 µg/ml	1.0 mg/mL	1.0 mg/mL
CBD-A	1244-58-2	100 µg/ml			1,000 µg/ml				
CBD-V	24274-48-4								
CBDV-A	31932-13-5								
CBG	25654-31-3								
CBG-A	25555-57-1								
CBN	521-35-7					100 µg/ml	1,000 µg/ml	1.0 mg/mL	1.0 mg/mL
CBN-A	2808-39-1								
THC / (-)-Δ9-THC	1972-08-3		3 µg/ml	1,000 µg/ml		100 µg/ml	1,000 µg/ml	1.0 mg/mL	1.0 mg/mL
THC-A	23978-85-0	100 µg/ml			1,000 µg/ml				
THC-V	31262-37-0								
THCVA	39986-26-0								
D8-THC	5957-75-5								
CBC	20675-51-8								
CBC-A	185505-15-1								
CBC-V	24274-48-4								
CBL	21366-63-2								
CBC-O	55824-09-4								
CBL-A	40524-99-0								
CBT	31508-71-1								

Cannabinoid Potency Standards

Multi Component CRMs:

Cannabinoid	CAS Number	3 Part Mix		4 Part Mix			5 Part Mix				
		LGC	Restek	Absolute Standards	Absolute Standards	Cayman	Absolute Standards	Absolute Standards	Cayman	LGC	LGC
		LGC-DRE-A50000201AL	RES-34014	AES-97896	AES-99361B	CAY-25106-1EA	AES-P98724K	AES-64382	CAY-25076-1EA	LGC-DRE-A50000189AL	LGC-DRE-A50000188AL
CBD	13956-29-1	1,000 µg/ml	1,000 µg/ml	100 µg/ml	100 µg/ml	1,000 µg/ml	100 µg/ml		500 µg/ml	1,000 µg/ml	500 µg/ml
CBD-A	1244-58-2						100 µg/ml		500 µg/ml	1,000 µg/ml	500 µg/ml
CBD-V	24274-48-4							100 µg/ml			
CBDV-A	31932-13-5										
CBG	25654-31-3				100 µg/ml						
CBG-A	25555-57-1							100 µg/ml			
CBN	521-35-7	1,000 µg/ml	1,000 µg/ml	100 µg/ml	100 µg/ml		100 µg/ml		500 µg/ml	1,000 µg/ml	500 µg/ml
CBN-A	2808-39-1										
THC / (-)-Δ9-THC	1972-08-3	1,000 µg/ml	1,000 µg/ml	100 µg/ml	100 µg/ml	1,000 µg/ml	100 µg/ml		500 µg/ml	1,000 µg/ml	500 µg/ml
THC-A	23978-85-0			100 µg/ml		1,000 µg/ml	100 µg/ml		500 µg/ml	1,000 µg/ml	500 µg/ml
THC-V	31262-37-0							100 µg/ml			
THCVA	39986-26-0										
D8-THC	5957-75-5					1,000 µg/ml		100 µg/ml			
CBC	20675-51-8							100 µg/ml			
CBC-A	185505-15-1										
CBC-V	24274-48-4										
CBL	21366-63-2										
CBC-O	55824-09-4										
CBL-A	40524-99-0										
CBT	31508-71-1										

Cannabinoid Potency Standards

Multi Component CRMs:

		5 Part Mix						6 Part Mix		
		LGC	LGC	LGC	LGC	NSI		Absolute Standards	Cayman	Cayman
Cannabinoid	CAS Number	LGC-DRE-A50000187AL	LGC-DRE-A50000186AL	LGC-DRE-A50000185AL	LGC-DRE-A50000184AL	NSI-C-928	NSI-C-929	AES-99361K	CAY-25077-1EA	CAY-35470-.5mL
CBD	13956-29-1	100 µg/ml						100 µg/ml	500 µg/ml	
CBD-A	1244-58-2	100 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml			100 µg/ml	500 µg/ml	1000 µg/ml
CBD-V	24274-48-4					1,000 µg/ml				
CBDV-A	31932-13-5		250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml				1000 µg/ml
CBG	25654-31-3							100 µg/ml	500 µg/ml	
CBG-A	25555-57-1		250 µg/ml	500 µg/ml	1,000 µg/ml		1,000 µg/ml			1000 µg/ml
CBN	521-35-7	100 µg/ml						100 µg/ml	500 µg/ml	
CBN-A	2808-39-1						1,000 µg/ml			
THC / (-)-Δ9-THC	1972-08-3	100 µg/ml						100 µg/ml	500 µg/ml	
THC-A	23978-85-0	100 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml			100 µg/ml	500 µg/ml	1000 µg/ml
THC-V	31262-37-0						1,000 µg/ml			
THCVA	39986-26-0		250 µg/ml	500 µg/ml	1,000 µg/ml		1,000 µg/ml			1000 µg/ml
D8-THC	5957-75-5						1,000 µg/ml			
CBC	20675-51-8					1,000 µg/ml				
CBC-A	185505-15-1					1,000 µg/ml				1000 µg/ml
CBC-V	24274-48-4									
CBL	21366-63-2					1,000 µg/ml				
CBC-O	55824-09-4									
CBL-A	40524-99-0									
CBT	31508-71-1									

Cannabinoid Potency Standards

Multi Component CRMs:

Cannabinoid	CAS Number	6 Part Mix				7 Part Mix			8 Part Mix		
		Cerilliant	Cerilliant	NSI	NSI	LGC	LGC	LGC	Cayman	Cerilliant	Cerilliant
		C-218-1ML	C-218S-1mL	C-925	C-924	DRE-A50000194AL	DRE-A50000195AL	DRE-A50000196AL	CAY-34043-0.5ML	C-219-1ML	C-219S-1mL
CBD	13956-29-1			1000 µg/ml	500 µg/ml				1000 µg/ml	500 µg/ml	500 µg/ml
CBD-A	1244-58-2	500 µg/ml	500 µg/ml	1000 µg/ml	500 µg/ml	250 µg/ml	500 µg/ml	1000 µg/ml			
CBD-V	24274-48-4								1000 µg/ml	500 µg/ml	500 µg/ml
CBDV-A	31932-13-5	500 µg/ml	500 µg/ml			250 µg/ml	500 µg/ml	1000 µg/ml			
CBG	25654-31-3			1000 µg/ml	500 µg/ml				1000 µg/ml	500 µg/ml	500 µg/ml
CBG-A	25555-57-1	500 µg/ml	500 µg/ml			250 µg/ml	500 µg/ml	1000 µg/ml			
CBN	521-35-7			1000 µg/ml	500 µg/ml				1000 µg/ml	500 µg/ml	500 µg/ml
CBN-A	2808-39-1										
THC / (-)-Δ9-THC	1972-08-3			1000 µg/ml	500 µg/ml				1000 µg/ml	500 µg/ml	500 µg/ml
THC-A	23978-85-0	500 µg/ml	500 µg/ml	1000 µg/ml	500 µg/ml	250 µg/ml	500 µg/ml	1000 µg/ml			
THC-V	31262-37-0								1000 µg/ml	500 µg/ml	500 µg/ml
THCVA	39986-26-0	500 µg/ml	500 µg/ml			250 µg/ml	500 µg/ml	1000 µg/ml			
D8-THC	5957-75-5								1000 µg/ml	500 µg/ml	500 µg/ml
CBC	20675-51-8								1000 µg/ml	500 µg/ml	500 µg/ml
CBC-A	185505-15-1	500 µg/ml	500 µg/ml			250 µg/ml	500 µg/ml	1000 µg/ml			
CBC-V	24274-48-4										
CBL	21366-63-2										
CBC-O	55824-09-4										
CBL-A	40524-99-0					250 µg/ml	500 µg/ml	1000 µg/ml			
CBT	31508-71-1										

Cannabinoid Potency Standards

Multi Component CRMs:

		9 Part Mix				10 Part mix				
		LGC	LGC	LGC	Restek	Cayman	LGC	LGC	LGC	LGC
Cannabinoid	CAS Number	DRE-A50000181AL	DRE-A50000182AL	DRE-A50000183AL	RES-34132	CAY-21305-1EA	DRE-A50000193AL	DRE-A50000192AL	DRE-A50000191AL	DRE-A50000190AL
CBD	13956-29-1	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBD-A	1244-58-2					250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBD-V	24274-48-4	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBDV-A	31932-13-5									
CBG	25654-31-3	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBG-A	25555-57-1					250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBN	521-35-7	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBN-A	2808-39-1									
THC / (-)-Δ ⁹ -THC	1972-08-3	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
THC-A	23978-85-0					250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
THC-V	31262-37-0	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml					
THCVA	39986-26-0						1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
D8-THC	5957-75-5	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml				
CBC	20675-51-8	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml	250 µg/ml	1,000 µg/ml	500 µg/ml	100 µg/ml	50 µg/ml
CBC-A	185505-15-1									
CBC-V	24274-48-4									
CBL	21366-63-2	1,000 µg/ml	500 µg/ml	250 µg/ml	1,000 µg/ml					
CBC-O	55824-09-4									
CBL-A	40524-99-0									
CBT	31508-71-1									

Cannabinoid Potency Standards

Multi Component CRMs:

Cannabinoid	CAS Number	11 Part Mix					14 Part Mix			16 Part Mix	
		Cayman	Cayman	Cayman	LGC	NSI	LGC	LGC	LGC	LGC	LGC
		CAY-32841-.5ML	CAY-21306-1EA	CAY-32842-1ML	DRE-A50000255AL	C-918	DRE-K50000202AL	DRE-K50000203AL	DRE-K50000204AL	DRE-K50000204AL	DRE-K50000207AL
CBD	13956-29-1	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBD-A	1244-58-2	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBD-V	24274-48-4	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBDV-A	31932-13-5						250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBG	25654-31-3	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBG-A	25555-57-1	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBN	521-35-7	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBN-A	2808-39-1										
THC / (-)-Δ9-THC	1972-08-3	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
THC-A	23978-85-0	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
THC-V	31262-37-0	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
THCVA	39986-26-0						250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
D8-THC	5957-75-5	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBC	20675-51-8	1 mg/ml	250 µg/ml	100 µg/ml	500 µg/ml	250 µg/ml	250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBC-A	185505-15-1										1,000 µg/ml
CBC-V	24274-48-4										
CBL	21366-63-2						250 µg/ml	500 µg/ml	1,000 µg/ml	1,000 µg/ml	1,000 µg/ml
CBC-O	55824-09-4										
CBL-A	40524-99-0										1,000 µg/ml
CBT	31508-71-1										

With over 100,000 products available, Jade Scientific offers an extensive selection to meet your needs. If you don't find what you're looking for in this catalog, please visit our website or contact a Jade Scientific sales representative for assistance.

New Lab Start-Up Checklist

CHEMICALS/REAGENTS:

- Acetone
- Acetonitrile
- Detergent
- Ethyl Acetate
- Hexane
- Media
- Methanol
- Nitric Acid
- Reagents
- Salts
- Silica Gel
- Sodium Chloride
- Solvents
- Standards
- Water

CHROMATOGRAPHY:

- Flash Media
- Syringes
- Vials

EQUIPMENT:

- Autoclaves
- Baths, Water
- Centrifuges
- Cryogenic Freezing
- Dispensers
- Freeze Dry Systems
- Homogenizers
- Hot Plates
- Incubators
- Ovens
- Forced Air
- General Purpose
- Vacuum
- Refrigerators/Freezers
- Rotary fs
- Shakers
- Benchtop
- Clamps
- Incubating
- Stirrers
- Thermal Cycler
- Ultrasonic Cleaners
- Vacuum Pumps
- Vortexers/Mixers
- Water Purification Systems

FURNITURE:

- Cabinets
- Carts
- Chairs/stools
- Countertops
- Drying Racks
- Faucets
- Fixtures
- Fume Hoods
- Glassware Washers
- Shelving
- Sinks
- Tables, Workstations

INSTRUMENTATION:

- Balances
- Barometers
- Conductivity
- Data Loggers
- Hydrometers
- Microplate Readers
- Microscopes
- Refractometers
- Spectrophotometer
- Timers
- Titration

LIFE SCIENCE:

- Cell Culture Supplies
- Microbiology
- Molecular Biology
- PCR
- Protein Purification

SAFETY:

- Apparel
- Chemical Storage
- Controlled Environment
- Ear Protection
- Emergency Response
- Eye Protection
- Eye Wash Kits
- Face and Head Protection
- First Aid
- Hand Hygiene
- Respirators
- Safety Identification

SUPPLIES:

- Aluminum Foil
- Apparel
- Aprons
- Caps
- Coveralls
- Cuff Tape
- Flocks
- Gowns
- Lab Coats
- Hoods & Face Covers
- Masks, Surgical
- Sleeves and Oversleeves
- Scrubs
- Shoe and Boot Covers
- Socks
- Autoclave Tape
- Bags
- Autoclave
- Biohazard Resealable Bags
- Bench Protectors
- Bottles
- Bunsen Burners
- Burets
- Centrifuge Tubes
- Clamps
- Clean Room Supplies
- Cleaners/Disinfectants
- Crucibles
- Face Masks
- Film
- Filtration
- Flasks
- Funnels
- Glassware
- Gloves
- Gas Regulators
- Hygiene
- Labeling Tape & Systems
- Liquid Nitrogen Dewars
- Microslides
- Parafilm
- Petri Dishes
- Pipets
- Pipettors
- Plasticware
- Plates
- Storage Boxes
- Thermometers
- Timers
- Tubes
- Vials
- Wash Bottles
- Wipe



Connect with us:

